

Compartment Pressure Measurement Traumatology




Compartment Pressure Measurement Basics

Definition

A compartment syndrome is a multi-factorial induced increase of the compartment pressure in closed and from fascies surrounded spaces, which lead to a disorder of the muscular microcirculation.

Due to the intracompartmental increase of the compartment pressure an acute or chronic blood flow disorder of the muscles occurs in association with an impairment of the muscle tissue and the related nerves.



Muscles are in closed, tight chambers - the muscle compartments - and each are surrounded from a fascie made of dense connective tissue.



Compartment Pressure Measurement Basics

Common affected Compartments


- Muscles of the lower leg especially in the area of the tibia
- Forearm muscles
- Muscles of the foot
- Muscles of the toes



Compartment Pressure Measurement Basics

Reasons for compartment syndrom

- Accumulation of fluids with haematoma and edema caused by bone fracture or contusion
- Permanent load of the muscles from long-distance runners
- Too fast muscle growth from strength athletes
- String impact with too tight plaster cast or high boots
- Long-lasting external pressure
- Drugs, medicines, alcohol
- Spontaneous bleedings with coagulopathy or administration of contrast media




Compartment Pressure Measurement Basics

Symptoms of compartment syndrom

- Feeling of tension*
- Neuronal deficits at a early stage*
- Muscular deficits*
- Strong pain*
- Palsy of the muscle group*
- Muscle stiffness
- Swelling of the muscle group

* Here is the feedback from the awake patient necessary!



Compartment Pressure Measurement Basics

Diagnosis of compartment syndrom

- Examination of neuronal and muscular deficits
- Imaging procedures provide especially at an early stage only unproductive statements.
- Measurement of intramuscular pressure IMP is very important with unconscious patients, because there is no feedback from the patients



Compartment Pressure Measurement Basics

Complications of the compartment syndrom

- Necrosis, mortify of muscle tissue
- Scarring of the muscels with a significant functio laesa right up to a ankylosis
- Damage of nerves leads to a neural deficit and permanent palsy





Compartment Pressure Measurement Basics

Conventional therapies

- Correct position of the patient
- Cooling of the affected muscle group
- Medication treatment
- Reduction of the wrong loading

Operative therapies

- Drainage for the diverting of fluid and with this for pressure relief
- Relief-OP in the form of a fasciotomy to improve the perfusion



Skin replacement

Compartment Pressure Measurement Basics

Prognosis of the process of a compartment syndrom

- Timely awareness and treatment are important
- Pressure relief within the first 4 hours
- Damages are mostly irreversible



Compartment Pressure Measurement Basics

Pressure limits


Compartment pressure = Intramuscular Pressure **IMP**

Simple pressure measurement (intramuscular pressure)

- < 10 mmHg Normal value for adults
- 10 – 20 mmHg Slight increased pressure
- 20 – 30 mmHg Impending compartment pressure syndrom
- > 30 – 45 mmHg Obvious compartment pressure syndrom

By determination of differential pressure Δp = diastolic blood pressure - intramuscular pressure

- ≥ 30 mmHg Normal pressure difference of adults
- < 30 mmHg Pathological pressure difference



Compartment Pressure Measurement Basics

Dependencies of the pressure measurement

- Position of the patient
- Position of the extremities
- The best measurement results can be achieved with comparison measurements in the healthy muscle




Figure 3. Anterior (A) and deep posterior (B) compartment IMPs as a function of knee and ankle position. (From Gershuni DH, Yaru NC, Hargens AR, et al: Ankle and knee position as a factor modifying intracompartmental pressure in the human leg. J Bone Joint Surg Am 66:1415–1420, 1984; with permission.)

Compartment Pressure Measurement

Measuring technique

Stryker Pressure Monitor:


Small pressure display device in combination with a fluid measurement via a cannula

Pro

- Continuous measurement is possible
- Economic: For each measurement only one new needle and a syringe filled with NaCl-solution is necessary
- Slot needle, thereby reduction of agglutination
- No sterilisation, because the cannula is disposable

Contra

- Injection of NaCl-solution leads to an increase of the pressure
- Needles agglutinate
- Air bubbles can falsify the measurement values
- No connection to a monitor possible
- Inconvenient handling



Compartment Pressure Measurement

Measuring technique

Stryker Pressure Monitor:

Small pressure display device in combination with a fluid measurement via a cannula

Pro

- Continuous measurement is possible
- Economic: For each measurement a new needle and a syringe filled with NaCl-solution is necessary
- Slot needle, thereby reduction of agglutination
- No sterilisation, because the cannula is disposable

Contra

- Injection of NaCl-solution leads to an increase of the pressure
- Needles agglutinate
- Air bubbles can falsify the measurement values
- No connection to a monitor possible
- Inconvenient handling



Compartment Pressure Measurement

Measuring technique

RAUMEDIC NPS3 and OMNIBAR E5F

Precise Measurement via a 5F microchip-catheter with piezo-resistive technique

Pro

- User friendly pressure display device NPS3 (or with NPS2 direct to the patient monitor or with MPR 1 DATALOGGER several measurements at the same time)
- Connection to the patient monitor with NPS2
- continuous measurement
- very easy handling
- no calibration of the catheter necessary
- catheter is a disposable product and steril supplied
- Connection of the measurement device (NPS3) with the implanted catheter is every time possible, that means one device (NPS3) can measure the pressure from different catheters alternating
- no hydrostatic artefacts
- no additional fluids which will be injected in the affected compartment



